

Adjusted by:

Tel:

Save the form to the desktop and open it with Adobe Reader.

E-mail:

Date:

The price changes depending on the choices made in the form **£:**

Standard configuration includes; Chassis, backrest, castors, rearwheel Spinergy 24", alu pushrim, Schwalbe One tire, one hand X brake, Carbon sideguards.

**Choice required:** Max. user weights : Seat width **33-45 cm 100 kg**,

**Optional** **Upcharge /Add-price** **Marked yellow** Click on "**blue underlined text**" for more info or image

**Original design**

**Alternative options**

**Other options**

**1) Panthera X**

Panthera X

**2) Chassis - Seat width** Art.no.

**33 cm      36 cm      39 cm      42 cm      45 cm**

Clear

**3) Backrest height 22-35 cm (**HH**)** Art.no.

**22 cm      23 cm      24 cm      25 cm      26 cm      27 cm      28 cm**  
**29 cm      30 cm      31 cm      32 cm      33 cm      34 cm      35 cm**

Clear

**4) Balance** Info

<b>Position 1</b> <b>Most</b> easy balanced	<b>Position 2</b> Less easy balanced	<b>Position 3</b> Medium balanced	<b>Position 4</b> Least easy balanced
--	---	--------------------------------------	--

The positioning of the backrest on the chassis affects the balance / COG of the wheelchair. "easy balanced" = easy to practice rear wheel balance. This can later easily be changed.

**5) Castor** Art.no.

**3" X** (87 mm) 90 mm S3

Clear

**6) Footrest**

Footrest Alu X  
36004BB

Footrest Titanium

36915BB

Carbon Footplate

65060BB

Footrest Titanium, extended. reinforced

36006BB needed for Freewheel

Plastic plate for footbow

36010BB

Clear

**7) Brake** Art.no.

One hand brake right One hand brake left

Clear

## 8) Rear wheel - without pushrim

[Rear wheel sizes](#)

[24" Spinergy X](#) black spokes

[25" Spinergy X](#) black spokes

[24" Spinergy X](#) white spokes

[25" Spinergy X](#) white spokes

Clear

[26" Spinergy X](#) black spokes

[26" Spinergy X](#) white spokes

## 9) Tyre - choose tyre

The price is for pair

[Schwalbe One](#) 24", 25", 26"

[Marathon Plus 24", 25", 26"](#)

[Rugged 24"](#)

[Right Run 24", 25", 26"](#)

[PU 24" 8 massive\)](#)

reTyre does not fit PU tyres

[reTyre skins 24"](#)

Clear

[reTyre skins 25"](#)

[reTyre skins 26"](#)

Clear

## 10) Pushrim

The price is for pair

[Aluminum X](#)

1904404 x2 (04=24", 05=25", 06=26")

Clear

[Titanium](#)

190DD11 (x 2

[Paragrip](#)

19051DDBT x 2

[Tetragrip](#)

19052DDBT x 2

[Maxgrip 22", 24"](#)

1903060 x 2 (60=24", 80=22"  
(contains latex)

DD= Wheeldiameter in inches

## Accessories

### Side guard Carbon [art.no.](#)

[Carbon without fender](#)

[Carbon small fender](#)

[Carbon large fender](#)

Clear

### Side guard S3 [video](#)

[Side guard S3](#)

### Cushion [Wedge cushion](#)

[5 cm](#)

572BB40

[Wedge cushion](#)

57710BB

## Arm hook

[Arm hook pair](#)

4200010

## Transport protection

[Housing](#)

9024009

[Rear end](#)

9024012

[Rear axle](#)

9024008

## Other

[Tube protection, front](#)

9900252 (2 pcs) [Tube](#)

[protection, rear](#)

9900252 (1 pcs)

[Freewheel](#)

Chose art.no. for Freewheel according to the footrest height of the wheelchair  
Read more about [Freewheel here](#)

## Set-up adjustments [Backrest angle](#)

Note measure (A) for desired back rest angle  
10-20 mm. Standard =15 mm

A= \_\_\_\_\_ mm

## [Footrest height](#)

Note measure (A) for desired footrest height

A= \_\_\_\_\_ mm

## Extra